



Insect specimen holder



rallekralle

[VIEW IN BROWSER](#)

updated 20. 3. 2022 | published 20. 3. 2022

Summary

To hold and smoothly rotate pinned bugs for microscopy or photography.

[Learning](#) > [Chemistry & Biology](#)

Tags: [stage](#) [bug](#) [microscope](#) [examination](#) [rotator](#)
[stereomicroscope](#)

Uses print-in-place springs to keep tension on the rotating parts.

Fill the hole in the wheel with some kind of clay and stick your pins in there.

You may have to adjust the size of the rotating cylinders until it's smooth. Grease might also be a good idea.

Model files

arm.stl





bug_rotator.stp



stand.stl



handle.stl



wheel.stl

License

This work is licensed under a
[Creative Commons \(International License\)](#)



Public Domain

- ✓ | Sharing without ATTRIBUTION
- ✓ | Remix Culture allowed
- ✓ | Commercial Use
- ✓ | Free Cultural Works
- ✓ | Meets Open Definition