



## Face Mask Assembly



Zack Slezak

[VIEW IN BROWSER](#)

updated 10. 4. 2020 | published 10. 4. 2020

### Summary

Print two of each piece. Put a piece of your preferred cloth between them and it will clamp them to hold the cloth.

[Healthcare](#) > [Home Medical Tools](#)

Tags: [healthandsafety](#) [covid19protection](#) [covidmask](#)

Print two of each piece. Put a piece of your preferred cloth between them and it will clamp them to hold the cloth. By using a shoe string or any other string you can create your own face mask. This is reusable and allows your to take the cloth out and wash it if need be.

### Print instructions

This model has been printed and tested with PET-G material.

### Model files



part-a.stl



part-b2.stl

## License ©

This work is licensed under a  
[Creative Commons \(International License\)](#)



**Public Domain**

---

- ✓ | Sharing without ATTRIBUTION
- ✓ | Remix Culture allowed
- ✓ | Commercial Use
- ✓ | Free Cultural Works
- ✓ | Meets Open Definition