



No Overhang Face Shield Bottom Reinforcement for Prusa Face Shield

 **Randompleb**

[VIEW IN BROWSER](#)

updated 3. 6. 2020 | published 3. 6. 2020

Summary

This remix fixes the overhang of the bottom reinforcement for Prusa face shield.

[Healthcare](#) > [Home Medical Tools](#)

Tags: [prusa](#) [face](#) [reinforcement](#) [faceshield](#) [ppe](#)
[covid19protection](#) [bottompiece](#) [prusafaceshield](#)
[covidfaceshield](#) [heathcare](#)

This remix fixes the overhang of the bottom reinforcement for Prusa face shield. Allowing you to print them at faster speeds without worrying about overhangs! It uses the same amount of filament as before but much easier to print and offers a little bit better bed adhesion.

Print instructions

Any layer height would work as there are no more overhangs. Tested with 200% speed with the standard 0.3mm PETG setting in Prusa slicer and it works perfectly fine. Your experience may vary.

This remix is based on



Prusa Face Shield

by Prusa Research

Model files



bottom-reinforcement-v2.stl

License

This work is licensed under a
[Creative Commons \(4.0 International License\)](#)



Attribution-NonCommercial

-
- ✗ | Sharing without ATTRIBUTION
 - ✓ | Remix Culture allowed
 - ✗ | Commercial Use
 - ✗ | Free Cultural Works
 - ✗ | Meets Open Definition