



Löwenstein prisma Filtergitter



Althron

[VIEW IN BROWSER](#)

updated 8. 12. 2023 | published 8. 12. 2023

Summary

Dieses Gitter verhindert das Herausfallen des Filters bei den prisma CPAP/APAP Geräten von Löwenstein.

[Healthcare](#) > [Home Medical Tools](#)

Tags: [cpap](#) [apap](#)

Durch Husten/Niesen kann der Filter an den Löwenstein prisma CPAP/APAP Geräten heraus gedrückt werden.

Dieses Gitter verhindert dies zuverlässig, ohne den Luftstrom zu behindern.


Gedruckt habe ich es aus PLA mit einer Layerhöhe von 0,12mm

Coughing/sneezing can cause the filter on Löwenstein prisma CPAP/APAP devices to be pushed out.

This grille reliably prevents this without obstructing the air flow.

I printed it from PLA with a layer height of 0.12mm

Model files



lowenstein_gitter.stl

License ©



This work is licensed under a
[Creative Commons \(4.0 International License\)](#)

Attribution-NonCommercial

-
- ✗ | Sharing without ATTRIBUTION
 - ✓ | Remix Culture allowed
 - ✗ | Commercial Use
 - ✗ | Free Cultural Works
 - ✗ | Meets Open Definition