



## ResMed CPAP P30i/ N30i Headgear Adjustment



ShadowDrakken

[VIEW IN BROWSER](#)

updated 3. 4. 2024 | published 3. 4. 2024

### Summary

CPAP mask sliding up off your head? This adjusts the headgear angle to cradle the base of the skull instead of the back.

[Healthcare](#) > [Home Medical Tools](#)

Tags: [mask](#) [frame](#) [cpap](#) [resmed](#) [headgear](#) [n30i](#) [p30i](#)

I got sick of my CPAP mask sliding up off my head every night because the headgear doesn't actually hold onto anything at the back of the head. This adjusts the angle downward so that the headgear can cradle the base of the skull and is unable to slide upward anymore.

You'll need to print two of them, with one mirrored, and they print best and have the best strength when printed on-edge with supports only to the print bed.

# Model files



resmed-cpap-p30i-n30i-headgear-adjustment-pair.stl



resmed-cpap-p30i-n30i-headgear-adjustment-single.stl

## License

This work is licensed under a  
[Creative Commons \(4.0 International License\)](#)



**Attribution—Noncommercial—Share Alike**

- 
- ✗ | Sharing without ATTRIBUTION
  - ✓ | Remix Culture allowed
  - ✗ | Commercial Use
  - ✗ | Free Cultural Works
  - ✗ | Meets Open Definition