



Bed calibration pattern for ANYCUBIC I3 MEGA



NewSpeak

[VIEW IN BROWSER](#)

updated 12. 4. 2022 | published 13. 2. 2022

Summary

Scroll test pattern can be makes Quick and Wide area test . 15 to 20 min in 25mm/s. pattern size : X260 Y190

[3D Printers](#) > [Test Models](#)

Tags: [anycubic](#) [anycubic4maxpro](#) [4maxpro](#)

Scroll test pattern can be makes Quick and Wide area test .

15 to 20 min in 25mm/s.

pattern size : X260 Y190

Category: 3D Printing Tests

Model files



[large_display_4maxprocalib_133955.stl](#)



4maxprocalib.stl

[Find source .stl files on Thingiverse.com](#)

License ©

This work is licensed under a
Creative Commons (4.0 International License)



Attribution

- ✗ | Sharing without ATTRIBUTION
- ✓ | Remix Culture allowed
- ✓ | Commercial Use
- ✓ | Free Cultural Works
- ✓ | Meets Open Definition