



# First layer calibration for Original Prusa i3 Mk3



JacobH4H3

[VIEW IN BROWSER](#)

updated 29. 3. 2022 | published 29. 3. 2022

## Summary

A great bed calibration for that oh so important first layer



1.76 hrs



1 pcs



0.15 mm



0.40 mm



PLA



11 g



Prusa  
MK3/S/S+

[3D Printers](#) > [Test Models](#)

A great bed calibration for that oh so important first layer. Intended for Prusa Original i3 MK3 but feel free to scale it to your bed size.

## Model files

**first-layer-calibration.stl**



# Print files

## first-layer-calibration.gcode



PLA 0.40 mm 0.15 mm 1.76 hrs 11 g Prusa MK3/S/S+

## License



This work is licensed under a  
[Creative Commons \(4.0 International License\)](#)

**Attribution—Noncommercial—Share Alike**

- 
- ✗ | Sharing without ATTRIBUTION
  - ✓ | Remix Culture allowed
  - ✗ | Commercial Use
  - ✗ | Free Cultural Works
  - ✗ | Meets Open Definition