



## Forearm crutch wall mount



Václav Blaschke

[VIEW IN BROWSER](#)

updated 22. 11. 2022 | published 22. 11. 2022

### Summary

Wall mount for one forearm crutch. PETG and a little more perimeters and infill are recommended.



4.88 hrs



1 pcs



0.20 mm



0.40 mm



PET



76 g



Prusa  
MK3/S/S+

[Healthcare](#) > [Medical Tools](#)

Tags: [crutchholder](#) [crutchesmount](#)

### Model files



[forearm\\_crutch\\_wallmount.stl](#)

# Print files



## forearm\_crutch\_wallmount\_02mm\_petg\_mk3s\_4h53m.gcode

PET 0.40 mm 0.20 mm 4.88 hrs 76 g Prusa MK3/S/S+

## License

This work is licensed under a  
[Creative Commons \(4.0 International License\)](#)



**Attribution-ShareAlike**

- 
- ✗ | Sharing without ATTRIBUTION
  - ✓ | Remix Culture allowed
  - ✓ | Commercial Use
  - ✓ | Free Cultural Works
  - ✓ | Meets Open Definition