



Test print

 **Adrian Črep**

[VIEW IN BROWSER](#)

updated 24. 4. 2023 | published 24. 4. 2023

Summary

Quick test print



0.25 hrs



1 pcs



0.20 mm



0.40 mm



PET



2 g



Ender 3 S1

[3D Printers](#) > [Test Models](#)

Tags: [print](#) [test](#) [testprint](#)

This is a short and quick test print to dial in your printer.

It prints in 15 minutes.

I used it to dial in retraction to minimize stringing and to improve smoothness on the vertical sections.

Model files



test.stl

Print files



test.gcode

⚙ PET ⚙ 0.40 mm ⚙ 0.20 mm ⌚ 0.25 hrs ⚖ 2 g

License ©

This work is licensed under a Standard Digital File Licence

Digital files have a strict non-commercial, personal use only license.

You shall not share, sub-license, sell, rent, host, transfer, or distribute in any way the digital file or 3D printed versions of this object, nor any other derivative work of this object in its digital or physical format. (including remixes of this object)

You can not host this files on other digital platforms, web stores or cloud repositories.

The objects may not be used without permission in any way whatsoever in which you charge money, collect fees.

-
- ✗ | No sharing or redistributing in any way of the 3D files or derivatives
 - ✗ | No remixing
 - ✗ | Non-commercial Use (only for personal use)
 - ✗ | Meets Open Definition